

Product datasheet for **TP726345**

Vinculin (VCL) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Vinculin/VCL
Species:	Human
Expression cDNA Clone or AA Sequence:	T331-L540
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Vinculin/VCL is produced by E. coli. The target gene encoding T331-L540 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7414
UniProt ID:	P18206
Synonyms:	Vinculin; VCL; Metavinculin (MV)
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Focal adhesion, Leukocyte transendothelial migration, Regulation of actin cytoskeleton



[View online »](#)