

Product datasheet for **TP726346**

UPF1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Regulator of nonsense transcripts 1/UPF1
Species:	Human
Expression cDNA Clone or AA Sequence:	L351-L654
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Regulator of nonsense transcripts 1/UPF1 is produced by E. coli. The target gene encodingL351-L654 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5976
UniProt ID:	Q92900
Synonyms:	Regulator of nonsense transcripts 1; ATP-dependent helicase RENT1; Nonsense mRNA reducing factor 1; Up-frameshift suppressor 1 homolog (hUpf1); UPF1; KIAA0221; RENT1
Protein Families:	Druggable Genome



[View online »](#)