

## Product datasheet for **TP726347**

### Munc 13 4 (UNC13D) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein unc-13 homolog D/UNC13D
Species:	Human
Expression cDNA Clone or AA Sequence:	A360-V540
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Protein unc-13 homolog D/UNC13Dis produced by E. coli. The target gene encodingA360-V540 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	201294
UniProt ID:	<a href="#">Q70J99</a>
Synonyms:	Protein unc-13 homolog D; UNC13D; Munc13-4



[View online »](#)