

Product datasheet for **TP726350**

Hexokinase 1 (HK1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Hexokinase-1/HK1
Species:	Human
Expression cDNA Clone or AA Sequence:	P76-A256
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Hexokinase-1/HK1 is produced by E. coli. The target gene encodingP76-A256 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3098
UniProt ID:	P19367



[View online »](#)