

Product datasheet for **TP726351**

CHERP Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human CHERP Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-L180
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human CHERP Protein is produced by E. coli. The target gene encoding A2-L180 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10523
UniProt ID:	Q8IWX8



[View online »](#)