

Product datasheet for **TP726361**

HSP90AA1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Heat shock protein HSP 90-alpha/HSP90AA1
Species:	Human
Expression cDNA Clone or AA Sequence:	E520-I755
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Heat shock protein HSP 90-alpha/HSP90AA1 is produced by E. coli. The target gene encodingE520-I755 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3320
UniProt ID:	<u>P07900</u>



[View online »](#)