

## Product datasheet for **TP726363**

### Kindlin 2 (FERMT2) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Fermitin family homolog 2/FERMT2
Species:	Human
Expression cDNA Clone or AA Sequence:	G260-L371
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Fermitin family homolog 2/FERMT2is produced by E. coli. The target gene encodingG260-L371 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10979
UniProt ID:	<a href="#">Q96AC1</a>
Synonyms:	Fermitin family homolog 2; FERMT2; Kindlin-2; Mitogen-inducible gene 2 protein (MIG-2); Pleckstrin homology domain-containing family C member 1 (PH domain-containing family C member 1); KIND2; MIG2; PLEKHC1



[View online »](#)