

## Product datasheet for **TP726366**

### **SAMHD1 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human SAMHD1 Protein
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	L320-V430
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	Recombinant Human SAMHD1 Protein is produced by E. coli. The target gene encoding L320-V430 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	25939
<b>UniProt ID:</b>	<u><a href="#">Q9Y3Z3</a></u>



[View online »](#)