

Product datasheet for **TP726373**

Cytokeratin 7 (KRT7) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Keratin, type II cytoskeletal 7/KRT7
Species:	Human
Expression cDNA Clone or AA Sequence:	V2-E710
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Keratin, type II cytoskeletal 7/KRT7is produced by E. coli. The target gene encodingV2-E710 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3855
UniProt ID:	P08729



[View online »](#)