

Product datasheet for **TP726374**

Lipoamide Dehydrogenase (DLD) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Dihydrolipoyl dehydrogenase, mitochondrial/DLD
Species:	Human
Expression cDNA Clone or AA Sequence:	F84-D205
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Dihydrolipoyl dehydrogenase, mitochondrial/DLD is produced by E. coli. The target gene encoding F84-D205 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1738
UniProt ID:	P09622



[View online »](#)