

Product datasheet for **TP726379**

Arp2 (ACTR2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Actin-related protein 2/ACTR2
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-E125
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol, 300mM NaCl, pH7.4
Note:	Recombinant Human Actin-related protein 2/ACTR2is produced by E. coli. The target gene encodingD2-E125 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10097
UniProt ID:	P61160
Synonyms:	Actin-related protein 2; Actin-like protein 2; ACTR2; ARP2
Protein Families:	Druggable Genome



[View online »](#)