

Product datasheet for **TP726380**

Glutamine Synthetase (GLUL) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Glutamine synthetase/GLUL
Species:	Human
Expression cDNA Clone or AA Sequence:	K241-C346
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Glutamine synthetase/GLUL is produced by E. coli. The target gene encodingK241-C346 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	2752
UniProt ID:	P15104



[View online »](#)