

## Product datasheet for **TP726382**

### SEPTIN5 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Septin-5/SEPTIN5
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-Q369
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol,300mM NaCl, pH7.4
Note:	Recombinant Human Septin-5/SEPTIN5 is produced by E. coli. The target gene encodingS2-Q369 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5413
UniProt ID:	<a href="#">Q99719</a>
Synonyms:	Septin-5; SEPTIN5; Cell division control-related protein 1 (CDCrel-1); Peanut-like protein 1; PNUTL1; SEPT5
Protein Families:	Druggable Genome
Protein Pathways:	Parkinson's disease



[View online »](#)