

Product datasheet for **TP726383**

ERK2 (MAPK1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Mitogen-activated protein kinase 1/MAPK1
Species:	Human
Expression cDNA Clone or AA Sequence:	S248-S360
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol,300mM NaCl, pH7.4
Note:	Recombinant Human Mitogen-activated protein kinase 1/MAPK1 is produced by E. coli. The target gene encoding S248-S360 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5594
UniProt ID:	P28482



[View online »](#)