

Product datasheet for **TP726384**

FHL1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Four and a half LIM domains protein 1/FHL1
Species:	Human
Expression cDNA Clone or AA Sequence:	K34-V135
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol,300mM NaCl, pH7.4
Note:	Recombinant Human Four and a half LIM domains protein 1/FHL1 is produced by E. coli. The target gene encodingK34-V135 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	2273
UniProt ID:	<u>Q13642</u>



[View online »](#)