

Product datasheet for **TP726385**

RPLP0 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Large ribosomal subunit protein uL10/RPLP0
Species:	Human
Expression cDNA Clone or AA Sequence:	V89-D317
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol,300mM NaCl, pH7.4
Note:	Recombinant Human Large ribosomal subunit protein uL10/RPLP0is produced by E. coli. The target gene encodingV89-D317 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6175
UniProt ID:	<u>P05388</u>



[View online »](#)