

Product datasheet for **TP726386**

PPP6C Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PPP6C Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-E223
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human PPP6C Protein is produced by E. coli. The target gene encoding A2-E223 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5537
UniProt ID:	O00743
Synonyms:	Serine/threonine-protein phosphatase 6 catalytic subunit; PP6C; PPP6C; PPP6
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)