

Product datasheet for **TP726388**

TOLLIP Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Toll-interacting protein/TOLLIP
Species:	Human
Expression cDNA Clone or AA Sequence:	E87-P274
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol,300mM NaCl, pH7.4
Note:	Recombinant Human Toll-interacting protein/TOLLIPis produced by E. coli. The target gene encodingE87-P274 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	54472
UniProt ID:	Q9H0E2



[View online »](#)