

Product datasheet for **TP726389**

LASP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LIM and SH3 domain protein 1/LASP1
Species:	Human
Expression cDNA Clone or AA Sequence:	N2-V210
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol,300mM NaCl, pH7.4
Note:	Recombinant Human LIM and SH3 domain protein 1/LASP1 is produced by E. coli. The target gene encoding N2-V210 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3927
UniProt ID:	Q14847
Synonyms:	LIM and SH3 domain protein 1; LASP1; LASP-1; Metastatic lymph node gene 50 protein (MLN 50); MLN50
Protein Families:	Druggable Genome



[View online »](#)