

Product datasheet for **TP726390**

PSME3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Proteasome activator complex subunit 3/PSME3
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-Y254
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol,300mM NaCl, pH7.4
Note:	Recombinant Human Proteasome activator complex subunit 3/PSME3is produced by E. coli. The target gene encodingA2-Y254 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10197
UniProt ID:	P61289



[View online »](#)