

Product datasheet for **TP726391**

RAB7 (RAB7A) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ras-related protein Rab-7a/RAB7A
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-S204
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol, 300mM NaCl, pH7.4
Note:	Recombinant Human Ras-related protein Rab-7a/RAB7A is produced by E. coli. The target gene encoding T2-S204 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7879
UniProt ID:	P51149



[View online »](#)