

Product datasheet for **TP726392**

SM22 alpha (TAGLN) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Transgelin/TAGLN
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-S201
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS,300mM NaCl,10%Glycerol, pH7.4
Note:	Recombinant Human Transgelin/TAGLN is produced by E. coli. The target gene encoding A2-S201 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6876
UniProt ID:	Q01995



[View online »](#)