

Product datasheet for **TP726393**

TMS1 (PYCARD) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PYCARD Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	L12-L162
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS,300mM NaCl,10%Glycerol, pH7.4
Note:	Recombinant Human PYCARD Protein is produced by E. coli. The target gene encoding L12-L162 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	29108
UniProt ID:	<u>Q9ULZ3</u>



[View online »](#)