

Product datasheet for **TP726396**

S100 alpha 6 (S100A6) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein S100-A6/S100A6
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-G90
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol,pH7.4
Note:	Recombinant Human Protein S100-A6/S100A6is produced by E. coli. The target gene encodingA2-G90 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6277
UniProt ID:	<u>P06703</u>



[View online »](#)