

## Product datasheet for **TP726398**

### Synaptopodin 2 (SYNPO2) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Synaptopodin-2/SYNPO2
Species:	Human
Expression cDNA Clone or AA Sequence:	V140-A380
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS,300mM NaCl,10%Glycerol, pH7.4
Note:	Recombinant Human Synaptopodin-2/SYNPO2is produced by E. coli. The target gene encodingV140-A380 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	171024
UniProt ID:	<a href="#">Q9UMS6</a>
Synonyms:	Synaptopodin-2; SYNPO2; Genethonin-2; Myopodin



[View online »](#)