

## Product datasheet for **TP726402**

### JMJD7-PLA2G4B Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cytosolic phospholipase A2 beta/PLA2G4B
Species:	Human
Expression cDNA Clone or AA Sequence:	Q659-H781
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS,300mM NaCl,10%Glycerol, pH7.4
Note:	Recombinant Human Cytosolic phospholipase A2 beta/PLA2G4Bis produced by E. coli. The target gene encodingQ659-H781 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	8681
UniProt ID:	<a href="#">P0C869</a>

[View online »](#)