

Product datasheet for **TP726406**

AP2M1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human AP-2 complex subunit mu/AP2M1
Species:	Human
Expression cDNA Clone or AA Sequence:	I2-I108
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human AP-2 complex subunit mu/AP2M1 is produced by E. coli. The target gene encoding I2-I108 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1173
UniProt ID:	<u>Q96CW1</u>



[View online »](#)