

## Product datasheet for **TP726407**

### Arp3 (ACTR3) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Actin-related protein 3/ACTR3
Species:	Human
Expression cDNA Clone or AA Sequence:	E39-E212
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol, 300mM NaCl,pH7.4
Note:	Recombinant Human Actin-related protein 3/ACTR3is produced by E. coli. The target gene encodingE39-E212 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10096
UniProt ID:	<a href="#">P61158</a>
Synonyms:	Actin-related protein 3; Actin-like protein 3; ACTR3; ARP3
Protein Families:	Druggable Genome



[View online »](#)