

Product datasheet for **TP726409**

TSG101 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Tumor susceptibility gene 101 protein/TSG101
Species:	Human
Expression cDNA Clone or AA Sequence:	G214-Y390
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Tumor susceptibility gene 101 protein/TSG101 is produced by E. coli. The target gene encoding G214-Y390 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7251
UniProt ID:	Q99816



[View online »](#)