

Product datasheet for **TP726410**

PRKACA Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PRKACA Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	D162-F351
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human PRKACA Protein is produced by E. coli. The target gene encoding D162-F351 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5566
UniProt ID:	P17612



[View online »](#)