

Product datasheet for **TP726412**

BUB3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Mitotic checkpoint protein BUB3
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-Y141
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol, 300mM NaCl,pH7.4
Note:	Recombinant Human Mitotic checkpoint protein BUB3is produced by E. coli. The target gene encodingT2-Y141 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	9184
UniProt ID:	<u>O43684</u>



[View online »](#)