

Product datasheet for **TP726413**

RACK1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Small ribosomal subunit protein RACK1
Species:	Human
Expression cDNA Clone or AA Sequence:	I31-W219
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol, 300mM NaCl,pH7.4
Note:	Recombinant Human Small ribosomal subunit protein RACK1 is produced by E. coli. The target gene encoding I31-W219 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10399
UniProt ID:	P63244



[View online »](#)