

Product datasheet for **TP726415**

GOLPH3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Golgi phosphoprotein 3/GOLPH3
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-L237
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol, 300mM NaCl,pH7.4
Note:	Recombinant Human Golgi phosphoprotein 3/GOLPH3is produced by E. coli. The target gene encodingT2-L237 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	64083
UniProt ID:	Q9H4A6
Synonyms:	Golgi phosphoprotein 3; GOLPH3; Coat protein GPP34; Mitochondrial DNA absence factor (MIDAS); GPP34



[View online »](#)