

Product datasheet for **TP726416**

STOM Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Stomatin/STOM
Species:	Human
Expression cDNA Clone or AA Sequence:	K55-V268
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol, 300mM NaCl,pH7.4
Note:	Recombinant Human Stomatin/STOM is produced by E. coli. The target gene encodingK55-V268 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2040
UniProt ID:	P27105
Synonyms:	Stomatin; STOM; Erythrocyte band 7 integral membrane protein; Erythrocyte membrane protein band 7.2; Protein 7.2b; BND7; EPB72
Protein Families:	Transmembrane



[View online »](#)