

## Product datasheet for **TP726422**

### GGPS1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Geranylgeranyl pyrophosphate synthase/GGPS1
Species:	Human
Expression cDNA Clone or AA Sequence:	K170-E300
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Geranylgeranyl pyrophosphate synthase/GGPS1 is produced by E. coli. The target gene encoding K170-E300 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	9453
UniProt ID:	<a href="#">O95749</a>



[View online »](#)