

Product datasheet for **TP726423**

Desmoplakin (DSP) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Desmoplakin/DSP
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-Y350
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Desmoplakin/DSP is produced by E. coli. The target gene encoding S2-Y350 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1832
UniProt ID:	P15924



[View online »](#)