

Product datasheet for **TP726426**

HINT2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Adenosine 5'-monophosphoramidase HINT2
Species:	Human
Expression cDNA Clone or AA Sequence:	G26-G163
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Adenosine 5'-monophosphoramidase HINT2is produced by E. coli. The target gene encodingG26-G163 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	84681
UniProt ID:	<u>Q9BX68</u>



[View online »](#)