

## Product datasheet for **TP726429**

### **PGAM1 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Phosphoglycerate mutase 1/PGAM1
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	A2-K254
<b>Tag:</b>	C-FLAG
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	Recombinant Human Phosphoglycerate mutase 1/PGAM1 is produced by E. coli. The target gene encoding A2-K254 is expressed with a FLAG tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	5223
<b>UniProt ID:</b>	<a href="#">P18669</a>



[View online »](#)