

## Product datasheet for **TP726430**

### HSBP1L1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human HSBP1L1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-T74
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human HSBP1L1 Protein is produced by E. coli. The target gene encoding D2-T74 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	440498
UniProt ID:	<a href="#">C9JCN9</a>



[View online »](#)