

Product datasheet for **TP726432**

CAPZB Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human F-actin-capping protein subunit beta/CAPZB
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-L243
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human F-actin-capping protein subunit beta/CAPZB is produced by E. coli. The target gene encodingS2-L243 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	832
UniProt ID:	P47756



[View online »](#)