

Product datasheet for **TP726434**

PABPC4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Polyadenylate-binding protein 4/PABPC4
Species:	Human
Expression cDNA Clone or AA Sequence:	P420-S644
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Polyadenylate-binding protein 4/PABPC4is produced by E. coli. The target gene encodingP420-S644 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	8761
UniProt ID:	<u>Q13310</u>



[View online »](#)