

Product datasheet for **TP726435**

RAB1A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ras-related protein Rab-1A/RAB1A
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-C205
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ras-related protein Rab-1A/RAB1A is produced by E. coli. The target gene encoding S2-C205 is expressed with a FLAG tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5861
UniProt ID:	<u>P62820</u>



[View online »](#)