

## Product datasheet for **TP726435**

### **RAB1A Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Ras-related protein Rab-1A/RAB1A
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	S2-C205
<b>Tag:</b>	C-FLAG
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	Recombinant Human Ras-related protein Rab-1A/RAB1A is produced by E. coli. The target gene encoding S2-C205 is expressed with a FLAG tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	5861
<b>UniProt ID:</b>	<a href="#">P62820</a>



[View online »](#)