

## Product datasheet for **TP726436**

### ARF4 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ADP-ribosylation factor 4/ARF4
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-N150
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human ADP-ribosylation factor 4/ARF4is produced by E. coli. The target gene encodingG2-N150 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	378
UniProt ID:	<a href="#">P18085</a>



[View online »](#)