

Product datasheet for **TP726437**

HMGB4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human High mobility group protein B4/HMGB4
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-S186
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human High mobility group protein B4/HMGB4is produced by E. coli. The target gene encodingG2-S186 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	127540
UniProt ID:	Q8WW32



[View online »](#)