

Product datasheet for **TP726438**

CALML3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Calmodulin-like protein 3/CALML3
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-K149
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Calmodulin-like protein 3/CALML3is produced by E. coli. The target gene encodingA2-K149 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	810
UniProt ID:	<u>P27482</u>



[View online »](#)