

Product datasheet for **TP726439**

Coatomer subunit delta (ARCN1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Coatomer subunit delta/ARCN1
Species:	Human
Expression cDNA Clone or AA Sequence:	K164-L511
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Coatomer subunit delta/ARCN1 is produced by E. coli. The target gene encoding K164-L511 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	372
UniProt ID:	P48444



[View online »](#)