

## Product datasheet for **TP726440**

### ADH5 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Alcohol dehydrogenase class-3/ADH5
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-C100
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Alcohol dehydrogenase class-3/ADH5 is produced by E. coli. The target gene encoding A2-C100 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	128
UniProt ID:	<a href="#">P11766</a>



[View online »](#)