

Product datasheet for **TP726441**

LRPPRC Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LRPPRC Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	L611-A907
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human LRPPRC Protein is produced by E. coli. The target gene encoding L611-A907 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10128
UniProt ID:	P42704



[View online »](#)