

## Product datasheet for **TP726442**

### hnRNP Q (SYNCRIP) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Heterogeneous nuclear ribonucleoprotein Q/SYNCRIP
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-T350
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Heterogeneous nuclear ribonucleoprotein Q/SYNCRIP is produced by E. coli. The target gene encoding A2-T350 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10492
UniProt ID:	<a href="#">O60506</a>



[View online »](#)